

## Zimmer Z Wave Hand Piece *HP<sup>Mini</sup> PT* For Smaller Areas

Introducing the new Zimmer Hand Piece, *HP<sup>Mini</sup> PT*, specifically designed for smaller, harder to reach areas; works with all Z Wave<sup>Pro</sup> and enPuls units. An ergonomic and smaller-than-standard Z Wave hand piece, the *HP<sup>Mini</sup> PT* allows for more focused and refined treatment. It is comfortable for both patient and practitioner and allows for efficient energy transfer from hand piece to tissue structures as deep as 35 mm.



The *HP<sup>Mini</sup> PT* includes the smaller hand piece and two smaller, differently shaped applicator heads. Their sizes are 6 mm and 15 mm (versus the standard 39 mm). They can be quickly changed without the use of tools, allowing treatments to be adjusted as needed.

The professional price is \$2,995 and the *HP<sup>Mini</sup> PT* is backed by Zimmer's 2 million pulse warranty.

**Set consists of: (1) Z Wave Mini PT hand piece, (1) 6mm applicator head, (3) 15 mm applicator heads and (10) blue caps**  
Compatible only with Z Wave<sup>Pro</sup>, not earlier versions of Z Wave. Inquire about upgrading your earlier Z Wave version.