

Z Wave^Q

with ultra-quiet,
water-cooled
hand piece
technology



Zimmer
MedizinSysteme

the "Q" is for "Quiet"
softer sound - equally effective

Zimmer



outperform

Radial Pulse Technology Was Never This Quiet - and Felt This Good

The new Z Wave^Q system is equipped with an innovative water-cooled generator ensuring high efficiency and durability.

The new ergonomic hand piece with "shock absorber" offers the industry's quietest operation, the highest level of safety, and maximum comfort to the patient.

Simple to Use

The large touch screen 12" display technology allows for quick selection of preset and saved programs. Application specific programs can be saved easily. Changes in treatment parameters are just a few taps away.

Reliable and Convenient

- Gentle, painless, and quick treatments
- No disposables
- No down time
- Non-invasive treatment
- Easy to operate
- Consistently reproducible results
- Visible results
- Quietest operation

outrun

Worry Free

All of this with a super 2-year limited warranty on the system and 2,000,000 pulses on the hand piece.

outsmart



Mode of Action

- The Zimmer Z *Wave*^o transmits powerful sound waves into tissue. The energy then spreads out radially to enhance the treatment effect.
- The body responds to the Z *Wave*^o sound waves with increased metabolic activity in the form of lymphatic drainage and increased blood circulation.
- The end result is an improvement in the appearance of cellulite and other body contouring procedures.



F=ma: Soft Shot Technology

Zimmer electromagnetic technology uses a projectile that has a 5 times greater **mass** than pneumatic systems. This makes for equal **Force** using 5 *times less* acceleration. Less **acceleration** means a softer, more comfortable pulse on the surface of the skin. We call this **Soft Shot Technology**.



Zimmer MedizinSystems
3 Goodyear, Suite B
Irvine, CA, 92618
1 (800) 327-3576
info@zimmerusa.com
www.zimmerusa.com
shop.zimmerusa.com

Zimmer
AestheticDivision

Technical Data

Technology	Compressor-free ballistic Radial Pulse Technology system with electromagnetic generator as a projectile accelerator. 1-22 Hz, 60-185mJ
Programs	5 preset options plus treatment protocols
Control Unit	12" color touch screen for all software operations, knobs for energy and frequency settings. Hand piece operation via multi-directional foot switch
Memory / Update	SD card for application menu memory, custom programs (99 spaces), error log, and firmware updates
Voltage	110 volts, 60 Hz
Dimensions	L 21" / W 21.5" / H 55"
Weight	43 lbs. (without hand piece)
Warranty	2,000,000 pulses on hand piece 2-year limited warranty on system

outdo