



Zimmer
MedizinSystems

PF Toner^{Pro}

Pelvic Floor Therapy

Kegels without the exercise

Passive pelvic floor muscle training with electromagnetic field energy for the treatment of male and female urinary incontinence.

May also contribute to additional benefits in intimate health and wellness.

More than 25 million adults, men and women, in the USA experience temporary or chronic urinary incontinence.*

Symptoms may have an effect on quality of life.

This may include urine leakage when coughing, sneezing or exercising. Some can experience a frequent/urgent need to urinate, and hesitancy to engage in sexual intimacy.

Magnetic field therapy treatment for incontinence induces a current which causes a comfortable series of pelvic floor muscle exercises.

This strengthening of the muscles leads to the improvement of incontinence disorders.

** Source: Johns Hopkins Medicine*



Pictured with optional chair

Pictured with optional chair



Powerful Therapy with *PF Toner*^{Pro}

Reverse the clock on incontinence

Therapy

Choose from preset programs or customize and save your own. The treatment lasts from 20 to 30 minutes. During the treatment, the patient sits in a comfortable chair while the pelvic floor muscles are optimally stimulated. The unique design allows for static and dynamic approaches, facilitating treatment from various vectors.

Recommendations

Typically, one weekly treatment for 8 to 10 weeks is recommended, depending on the patient's condition and response. Some patients report an improvement in incontinence after the first session.



No Consumables



No Disposables



Sensible ROI



New Revenue Stream



Non-invasive



Painless



Patient Stays Fully Clothed



Quick Lunchtime Treatment



Pictured with optional chair



Quality Design

German-made and crafted with high-quality fabric upholstery.



Easy To Use

Precise energy control. Setup is fast and easy, and treatment times typically last between 20-30 minutes.



Optional Chair with Interchangeable Inlays

(Figure A & B)

Custom design that includes 3 different seat positions for optimal treatment efficacy.



Static and Dynamic Treatment

Applicator can be used with and without the optional chair to treat from different vectors.



Optional Cushion

(Figure C)

Custom cushion fits your choice of seating.



Predefined Protocols

(Figure D)

Predefined programs for the treatment of urinary incontinence, as well as general pelvic floor training.



Customize & Save Favorites

Easily customize treatments and store "Favorites" in your protocol settings.



Treatment May Allow for Minimal Monitoring

The practitioner initially establishes duration and settings to ensure a comfortable and effective treatment for each patient.



Hygiene

The optional chair cover and optional cushion are easily cleaned following standard medical device cleaning protocols.

Specifications

1- 150 Hz Frequency	1 Channel	136 lb Approx. Tower	33" H 33" W 43" L PFToner ^{Pro} Tower	30" H 30" W 53" L PFToner ^{Pro} Optional Chair
MODES Predefined Programs AND Programmable Favorites		173 lb Approx. Optional Chair		

Zimmer MedizinSystems is in no way affiliated with Zimmer Biomet™

www.ZimmerUSA.com



Ⓐ Pictured with optional chair




Ⓑ Interchangeable inlays for optional chair



(bottom)

Ⓒ Optional cushion pictured with a typical office chair



Ⓓ Use predefined programs or save custom "favorites"

Zimmer MedizinSystems
3 Goodyear, Suite B
Irvine, CA 92618
(800) 327-3576 (office)
info@ZimmerUSA.com
www.ZimmerUSA.com
shop.ZimmerUSA.com

Zimmer
MedizinSystems